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Telefax (760) 598-6193

AW 1.14610 REV C

	CAUTION: FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH CARE PRACTITIONER.
	WARNING: THIS DEVICE IS NOT INTENDED FOR USE BY CHILDREN UNDER 12, DIABETICS, INCAPACITATED INDIVIDUALS, AND THOSE WITH DECREASED SKIN SENSITIVITY OR POOR CIRCULATION.
	WARNING: CAREFULLY READ POLAR CARE CUB AND COOLING PAD USE INSTRUCTIONS AND WARNINGS PRIOR TO OPERATION.
	WARNING: IF A NOTICEABLE CHANGE IN SKIN APPEARANCE IN THE AREA OF THE COLD APPLICATION IS OBSERVED, SUCH AS BLISTERING, DISCOLORATION, OR INCREASED SWELLING FOR MORE THAN AN HOUR AFTER USE, DISCONTINUE USE OF THIS PRODUCT AND CONSULT YOUR PHYSICIAN.

POLAR CARE CUB

CONTENTS:

1 POLAR CARE CUB, a localized cold therapy hand pump and ice water container for use with BREG cooling pads (sold separately).

USE INSTRUCTIONS

The Polar Care CUB must be used with a BREG Polar Pad, CUB model.

1 Apply Polar Pad according to application instructions included with it.

2 Fill cooler with ice to top. Add water to within 1/2 inch of top of cooler and screw on lid tightly. Firmly insert tubing stem into hole in cooler lid.

3 Connect pad by pushing couplings together firmly until you hear a “click.”

4 Place cooler on a table or stand at or above the height at which hand bulb will be used. Squeeze hand bulb slowly until red flow indicator ball is clearly visible in center of window. Do not squeeze hand bulb harder than necessary to keep ball centered.

5 It takes approximately 10 squeezes of hand bulb to completely fill pad. Once pad is full, stop pumping and allow water in pad to cool affected area for 5-10 minutes. Then squeeze hand bulb 3 to 4 times every few minutes or as necessary to maintain the cold effect. If skin temperature is cooler than comfortably tolerated, squeeze bulb less frequently. If skin temperature is warmer than desired, squeeze bulb more frequently.

6 Replenish ice and water before ice has melted completely.

7 To end session, stop pumping and remove pad. Drain as much water from system as possible by elevating pad above cooler and rolling it up slowly from end to squeeze out all remaining water. Empty ice and water and dry all components.

TIPS

Some condensation on tubes and pad is normal, particularly in humid environments. Most times that leaking is suspected, condensation is the actual cause. To reduce condensation at pad, provide insulating layers of material over pad before the securing wrap is applied.

Keep the height of cooler at or above the height of hand bulb in use. It is also preferable to keep hand bulb and cooler above pad height if possible.

If it is difficult to squeeze hand bulb and move the red ball, then check to be sure that there are no folds or kinks in pad when it was applied.

A thin dressing layer only should be applied between pad and skin. Do not apply pad over multiple layers of dressing or affected area will not cool properly.

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