

The Caregivers' Gait Belt



The ErgoBelt™ is the only gait belt designed by caregivers for caregivers. The ErgoBelt™ assists caregivers when walking a patient. Constructed with a sturdy nylon material and four double-stitched reinforced handles, this one size fits all ergonomically inclined device allows caregivers to assist properly when walking with a patient and boosts confidence in the patient's ability to move.

General Usage

The ErgoBelt™ allows caregivers to assist, when needed, those patients who can walk on their own, however, may need extra support. The belt is to be used only when walking with a patient. It is not intended for, or is used for, patient lifting or transferring.

This belt is designed for walking with a patient. The ErgoBelt™ will give the caregiver something to hold on to in case the patient encounters trouble in walking or may need a little extra support. The belt can encourage patients to walk more often, knowing that a nurse can be there to assist without placing any strain on the patient's body or clothing.

Important Prior To Usage

This brochure is only a guide, as alternative methods may be used for patient handling. We strongly recommend that your facility administrator or risk manager review and approve all protocols used for safe patient handling. Prior to using any equipment, be sure to check for signs of wear to ensure proper functionality. When using the ErgoBelt™ for the first time, practice with a fellow caregiver until you have a level of comfort with the product.


Cleaning Instructions

The ErgoBelt™ can be cleaned with a disinfectant. It can also be hand or machine washed up to 80°C. Do not tumble dry or wring. Air dry only. Inspect frequently for signs of wear. Withdraw from use if any faults are apparent and replace immediately

Correct Fitting Of Buckle

For proper buckling, be sure to first check that the straps are not twisted. Slide the buckle clip (A) into the buckle lock (B) until securely fastened. To remove the belt, simply press the sides of the buckle clip (A) together and slide outward.

Both Parts Of The Buckle Are Slightly Curved And Should Only Be Fitted Together One Way



Buckle should secure without force and undo easily by squeezing the two side prongs in the middle of the buckle.

Visit us online at www.ergosafeusa.com to see our full line of ergonomically designed products.

Item No. 2020-0022



ErgoSafe Products, LLC • 142 Firwood Drive, St. Louis, MO 63119
314-968-9066 • www.ergosafeusa.com