



Instructions for Use








---

UNLOADER<sup>®</sup> HIP


		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	12
ES	Instrucciones para el uso	15
NO	Bruksanvisning	18
DA	Brugsanvisning	21
SV	Bruksanvisning	24
FI	Käyttöohjeet	27
NL	Gebruiksaanwijzing	30
PL	Instrukcja użytkowania	33
JA	取扱説明書	36
ZH	中文说明书	39
KO	사용 설명서	42



11

12

13

Medical Device

## PRODUCT OVERVIEW (Figure 1)

- a. Pulley System
- b. Rotation Control Strap
- c. Thigh Strap
- d. Belt Arm
- e. Front Closure
- f. Zipper
- g. Upper Alligator Clip
- h. Lower Alligator Clip
- i. SmartDosing® Dial
- j. Trochanter Pad

## INDICATIONS FOR USE

Conditions benefitting from hip repositioning and/or compression around the greater trochanters, this may include:

- Hip osteoarthritis (OA)
- Femoroacetabular impingement (FAI)
- Labrum tears

## CONTRAINDICATIONS

None known.

## GENERAL SAFETY INFORMATION

**CAUTION:** For patients with peripheral vascular disease, neuropathy, and sensitive skin, physician supervision is recommended. If redness appears, decrease usage of the brace. Increase slowly as skin adapts.

**WARNING:** In case of functional change or functional loss, the patient should stop using the device and contact a healthcare professional.

**WARNING:** If the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

This device is for single patient use.

## FITTING INSTRUCTIONS

**NOTE:** At initial evaluation, fold belt arms and straps instead of cutting. When appropriate length has been determined, cut as needed to obtain final fit.

### **Rotation Control Strap Placement (Figure 2)**

- Before donning the Unloader Hip, fasten Rotation Control Strap on the affected side by threading the Upper Alligator Clip through a slit in belt arm, and attaching it to button hook.
- SmartDosing Dial should be at "0" mark.
- Temporarily fasten lower strap end at the bottom of trousers, so that it doesn't tear at the Lycra during fitting.
- Have user don the brace.

### **Belt Arm Adjustment (Figure 3–4)**

- Detach Front Closures from Belt Arms.
- Fold/Cut Belt Arms vertically, so that they are aligned with mid-line.
- Re-attach Front Closures.

- Close Belt Arms so that they overlap tightly but comfortably around waist.

### **Rotation Control Strap Adjustment (Figure 5–7)**

- Detach Upper Alligator Clip. Fold/Cut strap so that the SmartDosing edge is positioned below the Belt Arm line (approximately 2 cm/0.79”).
- Detach Lower Alligator Clip, and fold/cut strap to obtain ideal tension. It should feel tight, but not uncomfortable, when wrapped around the thigh and attached at outer side of leg.
- With leg rotated slightly outwards, set SmartDosing dial at “5” mark. When rotating leg towards front and lifting it slightly, patient should experience agreeable tension in leg.
- Repeat Lower Alligator Clip adjustment if necessary. The objective is to obtain adequate tension with SmartDosing at the “5” mark.

### **Thigh Strap (Figure 8–9)**

- Position Thigh Strap so that its dot is positioned over the dot on the Rotation Control Strap. Thread the slim end through the opening, and fasten securely.

### **Pulley System Adjustment (Figure 10–12)**

- Verify length of Pulley System strings by pulling at handles, and attaching on Belt Arms.
- If the handles overlap, shorten strings by using the spooling system located inside the handles.
- Fasten Pulley System so that it sits tightly and comfortably on Belt Arms.

### **Trochanter Pad (optional) (Figure 13)**

Apply pad when peak pressure reduction on greater trochanter is needed, and/or to increase sealing effect on labrum.

- Position Trochanter Pad on the greater trochanter of the affected hip by attaching it under Belt Arm.

### **Final check**

Verify that the brace is sitting correctly on body, with all parts accurately positioned. Ask the patient to walk around, sit down, and stand up, to ensure proper fit. Adjust as necessary, then proceed to cut straps and belt arms to ideal length.

## **MAINTENANCE**

- Hand wash using mild soap and rinse thoroughly. For machine wash, use a wash bag, and a program for delicate wash. Do not use fabric softener.
- Avoid prolonged contact with salt and chlorinated water. If such contact occurs, rinse well with tap water and air dry.
- Avoid prolonged exposure to strong sunlight.
- Air dry (do not tumble dry).

## **REPORT A SERIOUS INCIDENT**

### **Important notice to users and/or patients established in Europe:**

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## **LIABILITY**

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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- a. Zugsystem
- b. Rotationskontrollgurt
- c. Oberschenkelgurt
- d. Orthesen-Gurt
- e. Vorderer Verschluss
- f. Reißverschluss
- g. Obere Alligatorklemme
- h. Untere Alligatorklemme
- i. SmartDosing® Einstellrad
- j. Trochanter-Polster

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