medi Butler

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

medi Butler

Instructions for use

Step 1: Stand the medi Butler on a firm surface with the broader side of the base facing you, and the opening in the semicircular frame (stirrup) facing away from you.

Step 2: Lay the compression stocking in the stirrup with the toe of the stocking facing forwards. Turn the top of the stocking inside out and over the frame (Fig. 1).

Step 3: Slide the stocking down the stirrup until the heel pocket appears at the edge of the stirrup (Fig. 2).

Step 4: Place the medi Butler on the floor, if it is not there already. Slip your foot, toes first, into the toe of the stocking. It is recommended that you lean or hold onto something firm for stabilization (e.g. table, wall, or chair) (Fig. 3).

Step 5: Now step into the medi Butler. Adjust the toe and heel pocket material if necessary (Fig. 4)

Step 6: By alternating pulling and releasing the two handles of the medi Butler, slide the stocking up your leg towards the knee (Fig. 5).

Step 7: When the stocking is just below the knee, tilt the medi Butler slightly backwards and remove it (Fig. 6).

Step 8: To finish applying your stocking, slide the compression stocking evenly up your leg. Ease the stocking up bit by bit without pulling on the top band (Fig. 7)

Tips:

- Before donning your medi stockings, make sure that your legs are dry and remove all jewelry.
- Avoid pressing your fingernails into the stockings.
- Make sure that you never pull on the stocking's top band.
- The medi Butler can be cleaned with standard cleaning agents or disinfectants.

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