

509 Hallux Valgus Strap 510 Hallux Valgus Night Support

DE Gebrauchsanweisung	4
EN Instructions for use	7
FR Instructions d'utilisation	11
IT Istruzioni per l'uso	15
ES Instrucciones de uso	19
NL Gebruiksaanwijzing	23
SV Bruksanvisning	27
DA Brugsanvisning	31
NO Bruksanvisning	34
FI Käyttöohje	38
PL Instrukcja użytkowania	42
RU Руководство по применению	45

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Material	PA, PES, EA, PUR, PCM
-----------------	-----------------------

INFORMATION

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INFORMATION

Last update: 2017-04-26

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.

► Please keep this document in a safe place.

These instructions for use provide important information on applying the foot supports 509 hallux valgus strap and 510 hallux valgus night support.

2 Intended use

2.1 Indications for use

The support is intended **exclusively** for fittings of the lower limbs and **exclusively** for contact with intact skin.

The support must be used according to the indication.

2.2 Indications

- Hallux valgus (malposition of the metatarsophalangeal joint of the big toe)

Indications must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

Not known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; pronounced varicose veins, especially with return flow impairment, lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support will be applied; sensory and circulatory disorders in the foot area, e.g. in case of diabetic neuropathy.

2.4 Effects

The support promotes the muscular correction of toe positioning and alleviates pain.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

NOTICE

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

NOTICE

Contact with salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Application

CAUTION

Incorrect or excessively tight application

Skin irritation due to overheating, local pressure points due to tight application to the body

- ▶ Ensure that the product is applied properly and fits correctly.
- ▶ When in doubt, do not continue to use the product if skin irritation occurs.

509 hallux valgus strap

- > Open all hook-and-loop closures on the support.
- > The patient is sitting on the front edge of a chair.
- 1) Pull the loop over the big toe (see fig. 1).
- 2) Pull the longer strap underneath the foot at an angle to the outside and wrap it around behind the heel (see fig. 2).
- 3) Pull the longer strap in the direction of the big toe and fasten it to the short strap with the hook-and-loop (see fig. 3).
- 4) **Optional:** Shorten the strap to the desired length (see fig. 4).

510 hallux valgus night support

- > Open all hook-and-loop closures on the support.
- > The patient is sitting on the front edge of a chair.
- 1) Pull the loop over the big toe (see fig. 5).
- 2) Pull the longer strap underneath the foot at an angle to the outside and wrap it around behind the heel (see fig. 6).
- 3) Wrap the longer strap once around the shin (see fig. 7).
- 4) Pull the longer strap in the direction of the big toe and attach it to the short strap with the hook-and-loop.
- 5) **Optional:** Shorten the strap to the desired length (see fig. 8).

4.2 Cleaning

- 1) Fasten all hook-and-loop closures.
- 2) Recommendation: Use a laundry bag or net.
- 3) Wash the support in warm water at **40 °C** with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

6.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

6.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH