Rolyan® Adapt-It Thermoplastic Pellets:


A703-60 and A703-61

A

B

C

Adapt-It Thermoplastic Pellets can be used to modify or build up hand grips as well as to adapt equipment and attach outriggers.

PRECAUTIONS

- After heating thermoplastic pellets, check their temperature before handling them.

CARE AND CLEANING

- The completed device will lose its shape in temperatures over 135°F (57°C). It should be kept away from sources of heat such as ovens, hot water, open flames, and sunny car windows.
- The thermoplastic device can be cleaned with soap and lukewarm water.

INSTRUCTIONS FOR USE

The softened pellets adhere to plastic and metal surfaces. To prevent adhering, cover surfaces

with petroleum jelly or mineral oil prior to working with the softened pellets.

Softened pellets may feel sticky to dry skin. To avoid stickiness, apply hand lotion to skin prior to working with softened pellets. Softened pellets will bond to and remove nail polish.

1. Soften the pellets in hot water, using one of the following methods:
 - To use a heat pan (do not use a Hydrocollator® unit), heat water to a temperature of 150° to 160°F (65–70°C). Pour a mound of pellets into the corner of the heat pan (A), filling no more than one-third of the capacity of the heat pan with pellets.
 - To use a heat-resistant container, fill it no more than one-third full with pellets and pour boiling water over them.

2. Using a disposable stirring device, stir the pellets. As they soften, they become transparent and form a ball around the stirring device (B).
3. When all the pellets are soft, lift the stirring device to remove the ball from the water, and blot it on a towel.
4. Pull the ball from the end of the stirring device and knead it with fingers. Roll it in the palms of the hands to squeeze out any remaining water.
5. Form it into the desired shape (C).
6. Allow it to cool and harden. To accelerate hardening, immerse in cold water or spray with Cold Spray (Product No. 7454).

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Patterson Medical, Ltd.