

## Walker Boots

McKesson offers a line of high-quality soft goods for orthopedic treatment and rehabilitation. Our products provide comfort and security for patients recovering from orthopedic surgeries, injuries or accidents. Soft and comfortable to wear, they are designed to achieve an anatomically correct fit. The line includes shoulder and knee immobilizers, arm slings, wrist and forearm splints, and more. McKesson provides a combination of quality, integrity and value in every product for orthopedic treatment.

### Features and Benefits

- The soft nylon/foam liner has an integrated pump design making the pneumatic liner simple to inflate and adjust
- Lightweight, sturdy design with anatomically shaped uprights that help to facilitate a better fit
- Integrated strapping to aid in full circumferential support and compression\
- Cushioned inner/outer soles help to absorb shock upon heel strike
- High and low top versions available
- **Walker Boot Air:** The soft nylon/foam liner has an integrated pump design making the pneumatic liner simple to inflate and adjust
- Not made with natural rubber latex



#### Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

#### Never Compromise on Care<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

## ORDERING INFORMATION

MFR #	DESCRIPTION	TYPE	SIZE	MEN SIZE	WOMEN SIZE	QTY
155-79-95502	Walker Boot, Low Top	Non-Air	X-Small	2 - 4 in	3.5 - 5.5 in	1 ea
155-79-95503	Walker Boot, Low Top	Non-Air	Small	4.5 - 7 in	6 - 8 in	1 ea
155-79-95505	Walker Boot, Low Top	Non-Air	Medium	7.5 - 10.5 in	8.5 - 11.5 in	1 ea
155-79-95507	Walker Boot, Low Top	Non-Air	Large	10.5 - 12.5 in	11.5 - 13 in	1 ea
155-79-95508	Walker Boot, Low Top	Non-Air	X-Large	12.5+ in	13.5+ in	1 ea
155-79-95492	Walker Boot, Standard	Non-Air	X-Small	2 - 4 in	3.5 - 5.5 in	1 ea
155-79-95493	Walker Boot, Standard	Non-Air	Small	4.5 - 7 in	6 - 8 in	1 ea
155-79-95495	Walker Boot, Standard	Non-Air	Medium	7.5 - 10.5 in	8.5 - 11.5 in	1 ea
155-79-95497	Walker Boot, Standard	Non-Air	Large	10.5 - 12.5 in	11.5 - 13 in	1 ea
155-79-95498	Walker Boot, Standard	Non-Air	X-Large	12.5+ in	13.5+ in	1 ea
155-79-95522	Walker Boot, Low Top	Air	X-Small	2 - 4 in	3.5 - 5.5 in	1 ea
155-79-95523	Walker Boot, Low Top	Air	Small	4.5 - 7 in	6 - 8 in	1 ea
155-79-95525	Walker Boot, Low Top	Air	Medium	7.5 - 10.5 in	8.5 - 11.5 in	1 ea
155-79-95527	Walker Boot, Low Top	Air	Large	10.5 - 12.5 in	11.5 - 13 in	1 ea
155-79-95528	Walker Boot, Low Top	Air	X-Large	12.5+ in	13.5+ in	1 ea
155-79-95512	Walker Boot, Standard	Air	X-Small	2 - 4 in	3.5 - 5.5 in	1 ea
155-79-95513	Walker Boot, Standard	Air	Small	4.5 - 7 in	6 - 8 in	1 ea
155-79-95515	Walker Boot, Standard	Air	Medium	7.5 - 10.5 in	8.5 - 11.5 in	1 ea
155-79-95517	Walker Boot, Standard	Air	Large	10.5 - 12.5 in	11.5 - 13 in	1 ea
155-79-95518	Walker Boot, Standard	Air	X-Large	12.5+ in	13.5+ in	1 ea

McKesson Medical-Surgical  
 9954 Mayland Drive  
 Richmond, VA 23233  
[mms.mckesson.com/mckesson-brands](http://mms.mckesson.com/mckesson-brands)

Post-Acute Care and Home Health  
 800.347.2456

Home Medical Equipment  
 888.822.8111

Primary Care  
 866.625.2679

McKesson Brands Support  
 800.777.4908