



3M™ Transpore™ Surgical Tape

Transparent, perforated plastic tape that tears easily, works reliably.



1
*Easy to tear,
even with gloves*



2
Eye dressing



3
Peripheral IV tubing



4
General Dressing

Features and benefits

- Easy, straight, bi-directional tear
- Allows clinician to tailor the tape for various dressings and devices
- Porous
- Easy to handle with gloves
- Good adhesion to skin and tubing
- Transparent
- Hypoallergenic and latex free
- Water-resistant

Suggested applications

- Securing tubing and devices (e.g., catheters, IV tubing)
- Anchoring bulky dressings
- When tape width must be customized



3M™ Transpore™ Surgical Tape



3M's Skin Health Program provides you with a single source for all your skin health needs including the maintenance of healthy skin, the protection of skin at risk, and the treatment of damaged skin.

To learn more about 3M's Skin Health Program and the family of products, visit us at www.3M.com. Contact your 3M Health Care Sales Representative or call the 3M Health Care Customer Helpline at 1-800-425-3030.



Application recommendations:

- 1) Always apply to clean dry surface.
- 2) To maximize adhesion, gently press or rub from center out to edges of tape.
- 3) To promote skin health, do not stretch or apply with tension over areas where tape will be adhering directly to skin.
- 4) Refer to "3M Trouble-free Taping" brochure for additional usage recommendations.



Ordering Information 3M™ Transpore™ Surgical Tape

Cat.No.	Size	Rolls/Box	Boxes/Case
1527-0	½ in. x 10 yd./1,25 cm x 9,14 m	24	10
1527-1	1 in. x 10 yd./2,25 cm x 9,14 m	12	10
1527-1	2 in. x 10 yd./5,0 cm x 9,14 m	6	10
1527-3	3 in. x 10 yd./7,5 cm x 9,14 m	4	10



Health Care

Corporate Office - Bangalore : Concorde Block, UB City
24, Vittal Mallya Road, Bangalore - 560 001.
Tel : 080-2223 1414 Fax : 080 - 2223 1444
www.3M.com/in

Help
line **3M**
Health Care **1-800-425-3030**
Toll Free