

**#5130 Bed Pull-Up**

**Instructions:**

Thread white webbing strap around the bedframe at the foot of the bed. With buckle facing towards you, feed the white webbing strap through the bottom opening with teeth. Pull until desired length and then feed through top opening above teeth. Wrap Velcro closure around buckle.

**Dimensions:**

Webbing strap: 85" long.

Ladder pull-up: 45" long.

# CLOSED

# 5130

The STRAP is 47 inches Long




Beige loops 22.5" around  
+ 1" wide


White strap 25" long  
x 3/4" wide

This part  
is a lead

TO OPEN LIFT  
UP ON BAR  
AND THEN THE  
2 TEETH WILL  
LIFT UP ALSO  
TO BE IN THE  
"OPEN" STATE  
SEE "OPEN PAGE



CLOSED



FRED SAMMONS, INC.

PART# SIZE

60761

5130 00

MADE IN U.S.A.

OPEN

← OPEN

TO OPEN  
LIFT ON BAR  
SEE "CLOSED"  
PAGE

FRED SAMMONS, INC.

PART#	SIZE
60761	33
5130	00

MADE IN U.S.A.