

Comfy™ Hand Contracture Management



Comfy GRIP

- Two rolls encourage finger extension of the MCP and PIP joints
- Great for contracted fingers
- Holds hand in functional C position
- Headliner Comfy cover

Comfy Finger Extender

- Frame designed to slide between palm and fingers of even the most contracted hands
- Two graduated rolls gradually increase extension of the MCP and PIP joints
- Ambidextrous



Comfy Air Hand Orthosis

- Gentle air bladder inflated as needed
- Gradual extension
- Side wings & Easy positioning
- Headliner Comfy cover
- Ambidextrous



Comfy Air Hand Roll

- Inflatable roll can easily slide through contracted hands by pulling on attached string
- Gentle air bladder inflated as needed
- Optional finger separator available
- Fits onto most Comfy Splints or use alone



Comfy Splints
We Make The World Comfy

info@lenjoymedical.com
www.comfysplints.com

13112 Crenshaw Blvd,
Gardena, CA 90249
T: 800-582-5332