

Eko DUO

Digital Stethoscope + One-Lead ECG.
For the complete picture.



Gain Unprecedented Insight Into Cardiac Function

The first device combining ECG and an electronic stethoscope, the DUO enables clinicians to get a comprehensive view of cardiovascular function without a costly, time-consuming work-up.



FDA
Cleared



Digital
stethoscope



1-lead
ECG



60x sound
amplification



Wireless
Bluetooth LE




Pairs with free
Eko App

Eko



DUO ECG + Digital Stethoscope

\$349 | DUO101
MSRP | MPN



Powerful Software Toolkit iOS, Android & Windows

Capture it


Save 15, 30, 60 or 120-second sound and ECG recordings for playback and visualization of waveforms.

Analyze it

Eko is developing advanced machine learning algorithms that analyze heart sounds and ECG tracings to help screen for heart disease.

Share it

Share recordings with colleagues for a second opinion or Livestream cardiology-grade Sounds and ECGs with an Eko enterprise plan.



*Eko's algorithms are in development and have not been cleared by the FDA.
TM and © 2019 Eko Devices, Inc. All Rights Reserved | LBL 070 Rev A