


# THE HANDCUFF™



## TIME TO GET HANDCUFFED?

HAND-BASED ORTHOSIS FOR THE TREATMENT OF METACARPAL FRACTURES

### 1010 HANDCUFF™

#### INDICATIONS:

- STABLE METACARPAL FRACTURES

#### FEATURES:

- 3 SIZES (S, M, L); FITS RIGHT OR LEFT
- CUSTOM MOLDED FOR INTIMATE FIT
- LOW PROFILE DESIGN LEAVES FINGERS AND WRIST FREE ENCOURAGING FUNCTIONAL USE
- STEP DOWN IMMOBILIZATION FOLLOWING TKO OR CAST REMOVAL
- PROMOTES FUNCTIONAL USE OF THE HAND WHILE PROTECTING HEALING FRACTURES



#### DESIGNED IN COLLABORATION WITH:



STEVEN SHIN, MD

#### ORDERING INFORMATION:

ORDER #	SIZE	MIN IN	MAX IN	MIN CM	MAX CM
1010	S	6.5	8.5	17	22
1010	M	8.5	9.5	22	24
1010	L	9.5	11	24	28

DIRECTIONS: MEASURE HAND CIRCUMFERENCE BELOW THE KNUCKLES

CALL TODAY! (800) 654-3241  
EMAIL: [ORDER@HEL-Y-WEBER.COM](mailto:ORDER@HEL-Y-WEBER.COM)



© 2018 HELY & WEBER



**HEL-Y-WEBER**  
BRACING INNOVATIONS

# THE HANDCUFF™ I-PLUS



**TIME TO GET HANDCUFFED?**

HAND-BASED ORTHOSIS FOR THE TREATMENT OF BOXERS FRACTURES



## 1048 HANDCUFF™ I-PLUS (INTRINSIC PLUS)

### INDICATIONS:

- STABLE METACARPAL FRACTURES
- ULNAR NERVE INJURIES (ANTI-CLAW ORTHOSIS)

### FEATURES:

- 2 SIZES (S/M, L/XL); RIGHT AND LEFT
- NEW COMPRESSION MOLDED BODY AND DORSAL STAY FOR INTIMATE FIT AND DESIRED POSITIONING
- STEP DOWN IMMOBILIZATION FOLLOWING TKO OR CAST REMOVAL
- DORSAL BLOCK PREVENTS FINGER EXTENSION WHILE ALLOWING FULL FINGER FLEXION
- OPTIONAL FINGER EXTENSION STAY WITH STRAP FOR CUSTOMIZABLE IMMOBILIZATION



### DESIGNED IN COLLABORATION WITH:



JAMES SAUCEDO, MD  
NOT PICTURED:  
ABEL CAVAZOS, OTR, CHT

### ORDERING INFORMATION:

ORDER #	SIZE	MIN IN	MAX IN	MIN CM	MAX CM
1048	S/M	6.5	8.5	17	22
1048	L/XL	8.5	10	22	25

SPECIFY LEFT OR RIGHT  
DIRECTIONS: MEASURE HAND CIRCUMFERENCE BELOW THE KNUCKLES

CALL TODAY! (800) 654-3241  
EMAIL: [ORDER@HELW-WEBER.COM](mailto:ORDER@HELW-WEBER.COM)



© 2018 HELY & WEBER

**HELW & WEBER**

BRACING INNOVATIONS