## **OXYGEN CYLINDERS**



## **ALUMINUM MEDICAL CYLINDERS**

- Brushed aluminum with clear coat protective finish
- 25-40% lighter than steel
- 5 year hydrostatic test
- White top cylinders available for international markets



















ALUM. SIZE	M4	M6	ML6	C [M9]	D	JUMBO D	E	M60	MM
CAPACITY (L)	106 L	161 L	165 L	244 L	398 L	623 L	670 L	1699 L	3455 L
CYL PER PALLET	144	144	81	81	81	72	81	42	25
WEIGHT (LBS)	1.7	2.3	2.8	3.6	5.1	8.4	7.5	22.3	40.7
CGA 870 STD POST VALVE	#110-0810	#110-0110	#110-0510	#110-0210	#110-0310	#110-0350	#110-0410	N/A	N/A
CGA 870 TOGGLE VALVE	#110-0820	#110-0120	#110-0520	#110-0220	#110-0320	#110-0351	#110-0420	N/A	N/A
CGA 540 VALVE & CARRY HANDLE	N/A	#110-0660	#110-0760						

 $\label{thm:cylinders} \textbf{Cylinders are available with "white" top for international markets - call for details. Weights and capacities are approximate.}$ 

## **HOME FILLING CYLINDERS**

- Genuine home filling valves and cylinders
- 100% tested, market proven design
- Avoid the hassle and expense of an integrated conserver and its repairs! Add a Respond C5 conserver to your home filling purchase.
- Home filling carry cases available, reference page 17





#111-3001

## **HOME FILLING CYLINDER VALVES**

#111-3000 Home filling valve

(Compare to HF2POST-R)

#111-3001 Home filling valve - D & E cylinders

(Compare to HF2POST-D)



ALUM. SIZE	ML6	C [M9]	D	E
CAPACITY (L)	165 L	244 L	398 L	670 L
CYL PER PALLET	81	81	81	81
CGA 870 HOME FILLING VALVE W/ TOGGLE	#111-0520	#111-0220	#111-0320	#111-0420
COMPARE TO:	HF2POST6	HF2POST9	HF2POSTD	N/A
•		•	•	