

# **DC365 Surgical N95 Respirator**

Surgical N95 Respirator **Technical Data Sheet** 

## **Product Description**

The Honeywell DC365 is a filtering facepiece Surgical N95 Respirator with nose clip. The DC365 surgical respirator does not contain latex or cellulose.

### Certification

This respirator is approved and cleared for use in healthcare settings under the FDA/NIOSH MOU 225-18-006. Standard Test Procedures used by NIOSH NPPTL under NIOSH Title 42 Code of Federal Regulations, Part 84. FDA registered. This Surgical N95 provides protection against particulates at the N95 filter efficiency level and conforms to recognized standards for flammability, fluid resistance and biocompatibility.

**TC Number** TC-84A-9250 DC365N95HC

Part Number

Description Surgical N95 Respirator

**Specifications** 

Fluid-resistance ASTM F1862 160 mmHg

Color White and green

**Description Material** 

Outer Laver Nonwoven filter media, dye-free and non-sensitizing Middle Layer Meltblown filter media Hydrophobic spunbond polypropylene Inner Laver Head Straps Polyester and spandex, ultrasonic-bonded Galvanized steel wire encased in polypropylene, ultrasonic-bonded Nose Clip

### Ordering Information

Part Number	Packaging	Bx,Bg/Case	Shipping Wt.	Dimensions
DC365N95HC	20/box	10 bxs/case	6.6lbs (3.0kg)	27x9.63x9.63 in (68.58x24.45x24.45 cm)

**Country of Origin** United States

**Technical Support** 

800-873-5242

See package labeling for additional warnings and precautions.