

# PHILIPS

OptiChamber Diamond

Valved holding chamber



## Designed to **improve medication delivery**

The anti-static material inside the OptiChamber Diamond allows the aerosol to be suspended for longer, giving your patients more time to inhale.<sup>1</sup>

Low resistance valves open freely, even under low paediatric flow rates, allowing the user to breathe easily through the chamber.

#1 pharmacy preferred brand of Spacer<sup>2</sup>



## Helping them stay in control



### Breathe easier

One-way valve stops the patient accidentally exhaling into the chamber.



### Detachable mask

The LiteTouch facemask provides a soft cushion that rests gently on the face.



### Device feedback

A flow alert 'whistle' tells patients when they are breathing in too quickly.

OptiChamber Diamond is designed for your patients, including those with a number of specific needs:



#### Helps coordination

The anti-static material suspends the aerosol longer, meaning your patient has more time to inhale their medication.<sup>1</sup>



#### Easy to maintain

The OptiChamber Diamond is easy to disassemble and can be easily cleaned with warm water and mild liquid detergent.



#### One size works for all

The stepped mouthpiece makes the OptiChamber suitable for use by children and adults without any configuration.



#### Supports caregiver

The exhalation valve visually confirms exhalation, helping the caregiver count patient breaths and coordinate actuation.



## Product specification table

Length	14.2 cm
Warranty	1 year
Volume	140 ml
Mouthpiece	Interfaces with 22 mm standard connectors
Materials	Chamber is anti-static ABS plastic, LiteTouch mask is polycarbonate (PC) silicone, valves are latex-free silicone

**#1** pharmacy preferred brand of Spacer.\*

\* Among retail pharmacies who claim to have a preferred brand of spacer in their pharmacy. Online survey conducted on 205 registered retail pharmacies.

1. Mitchell, J., et al. "Electrostatics and Inhaled Medications: Influence on delivery via pressurized metered-dose inhalers and add-on devices," *Resp. Care*. March 2007, Vol 52 No 3.

Find out more – [www.philips.com/respironics](http://www.philips.com/respironics)



© 2018 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Respironics Inc. and Respironics Respiratory Drug Delivery (UK) Ltd are Philips companies.



Respironics Respiratory Drug Delivery (UK) Ltd  
Chichester Business Park, City Fields Way  
Tangmere, Chichester PO20 2FT UK  
[www.philips.com/respironics](http://www.philips.com/respironics)

MCI 4108340 PN 1139760  
OptiChamber Diamond DataSheet, US Eng  
EB 04/07/18