

NEW!


PrimaFit™

External Urine Management System for Females

Address the #1 CAUTI risk factor, the indwelling catheter

PrimaFit is an external urine management system designed to help manage urinary incontinence in female patients and promote early catheter removal, addressing the #1 risk factor of catheter associated urinary tract infections (CAUTI). PrimaFit serves as an alternative to incontinence pads, diapers, or at times, catheters, while transitioning patients to independent continence.

When properly positioned, PrimaFit's ultra-soft wicking fabric absorbs and diverts urine away from the patient's skin. Urine is then absorbed into the system's core and suctioned into a collection canister, helping keep the skin dry and addressing a risk factor of incontinence-associated dermatitis (IAD).



No catheter. No CAUTI.™



Easy connection
Tiered Connection Port allows for easy connection to hospital tubing

Gentle, stays in place
Sure-Stay Silicone Adhesive Pad keeps device in place, easy and gentle removal

Helps keep skin dry
Ultra-soft wicking fabric absorbs and diverts urine away from skin

Flexible fit
Flex-Fit Core contours to patient anatomy and maintains shape for the duration of use

Stays in place
Tapered End Cap gently fits in perineal area to secure device in position

Approximately **75%** of hospital-acquired UTIs are associated with an indwelling urinary catheter.¹



PrimaFit
(1) PrimaFit
40 devices/case
Reorder #5400



Clean-up Cloths with PrimaFit
(1) package of 3 medium weight clean-up cloths
(1) PrimaFit
40 systems/case
Reorder #5452